

MEG II実験のための 背景ガンマ線同定用低運動量 陽電子タイミングカウンターの 実機製作と性能評価


Development status of low momentum positron timing counter
to identify BG gamma ray from radiative muon decay
in MEG II experiment

中浦正太、他MEG IIコラボレーション
東京大学
日本物理学会 2015年秋季大会
@大阪市立大学(杉本キャンパス)

- MEG II Experiment
- Radiative Decay Counter (RDC)
- Status & Schedule
- Scintillator Test
- MPPC Grouping
- Counter Production & Test
- Summary & Prospect

◇ $\mu \rightarrow e\gamma$ decay search

- In standard model with neutrino oscillation, the branching ratio is too small to detect the decay: $\sim O(10^{-54})$
- On the other hand, some beyond standard models predict large branching ratio: $\sim O(10^{-14} - 10^{-13})$ whose level our experiment can reach!!
- **Discovery of $\mu \rightarrow e\gamma$ = Discovery of new physics**



◇ **Signal Event** : 2-body decay from a muon at rest

- Both e^+ and γ have a **monochromatic energy** (52.8 MeV)
- They are emitted **time-coincidentally**, and
- **back to back**


Basic Idea

high precision measurement of energy, timing and opening angle of e^+ and γ

◇ Upgrade of MEG Experiment

T.Iwamoto(27aSN-8), S.Ogawa(27aSN-10), K.Ieki(27aSN-9)

Status:
Y.Uchiyama
(25aSG-2)



MEG I
Upper Limit:
 5.7×10^{-13}
(90% C.L.)

Phys. Rev. Lett. 110, 201801
- Published 13 May 2013



MEG II
Sensitivity
Goal:
 4×10^{-14}

S.Nakaura
(28aSG-7),
R.Iwai
(28aSG-8)


M.Nishimura(26aSN-1), M.Nakao(26aSN-2), K.Yoshida(26aSN-3)

- Doubled beam rate: 7×10^7 muons/s
- PMTs in LXe Detector are partially replaced with MPPCs
- A low mass stereo drift chamber
- A multi-tile scintillation timing counter
- New BG tagging detector is introduced

◇ Main Background Event : accidental background



- Radiative Decay Counter (RDC) tags Radiative Muon Decay actively (→Next page)
- 16 ~ 28% higher sensitivity can be achieved by introducing RDC





	sensitivity
without RDC downstream	5.0×10^{-14}
with RDC downstream	4.3×10^{-14}
with full RDC	★ $3.9 \sim 4.1 \times 10^{-14}$

16~28%


★ depending on detection efficiency (50%~80%)

◇ Active tagging of Radiative Muon Decay

- Radiative Muon Decay with a high energy (> 48 MeV) gamma ray **simultaneously emits** a low momentum positron (typically: ~ 5 MeV).




- It has small turning radius
→ Small counters are set on the beam axis
- measure **time coincidence** between γ and low momentum e^+
- measure **energy** of e^+ to distinguish Radiative Decay from Michel decay ← downstream only




◇ Upstream RDC

- Made of 704 plastic scintillation fibers
 - thickness = 250 μm
 - small effect on μ^+ beam transportation
- Separate μ^+ from e^+ using difference of energy deposit
- Measure **time coincidence** between γ and e^+ (from radiative decay)





The effect on the μ^+ beam properties ?

- In simulation,
 - stopping efficiency will be decreased by 3 %:
71.9% (w/o RDC) \rightarrow 69.8% (w/ RDC)
 - A small influence on the beam spread (6~7 %).
sigma X: 1.11 cm (w/o RDC) \rightarrow 1.19 cm (w/ RDC)
sigma Y: 1.07 cm (w/o RDC) \rightarrow 1.13 cm (w/ RDC)
- We plan to measure it with beam using a mock-up detector in this Autumn.




◇ Downstream RDC

- Already approved by the collaboration
- 12 Plastic scintillation bars
 - length: 7 ~ 19 cm, width: 1 ~ 2 cm, thickness: 5 mm
 - measure time coincidence between γ and e^+
- 76 LYSO crystals
 - $2 \times 2 \times 2 \text{ cm}^3$
 - measure energy deposit and distinguish Radiative Decay and Michel Decay




plastic scintillator



LYSO



This talk



◇ Plastic scintillator

Good performance within limited space


- BC-418 (Saint-Gobain)
 - large light yield
 - fast response (decay constant: 1.4 nsec)
- Two types of width: 1 cm & 2 cm
 - 1 cm: for reducing pileup events
- Series connected MPPCs
 - sharpen waveforms
- PCBs: compact readout circuits
 - MMCX cables are used.
 - 3 narrow bars are attached on a long PCB
- Required time resolution ~ 100 psec



◇ Plastic scintillator

setup




- 12 scintillation bars were tested.
- MPPCs were attached on both sides of a scintillation bar.
- 2 or 3 MPPCs were connected in series.
- optical grease was used for coupling.
- radioactive source: ^{90}Sr .
- trigger: hit in reference counter.



- calculate $[T(\text{ch1})+T(\text{ch2})]/2 - T(\text{reference})$
 - constant fraction time
- Fit the distribution with Gaussian function and the **sigma** is defined as time resolution



• Result: 120-145 psec

- This value includes resolution of the reference counter. (the resolution of reference counter : ~100 psec)
- The resolution of the bars is **good enough**.




◇ MPPCs for plastic scintillator

- Grouping
 - there are 12 scintillation counters
→ MPPCs are divided into 24 groups
 - measure I-V curve for all of MPPCs
 - MPPCs which show close I-V shapes are put in the same group.





I-V curves of 3 MPPCs in the same group

why?




◇ PCB design

- Consists of 4 layers:
 - shield layers on both sides for noise suppression



capacitor

- Capacitors for noise reduction
 - connecting ground of MPPCs with chassis ground (this is effective on reducing noise by our experience)
 - avoid ground loop.




◇ Assembling

- Conductive epoxy (CW2400)
 - attach MPPCs to a PCB without heating
 - strong mechanical bonds
 - good electrical conductivity




- Alignment jig




POM parts



- a π - shaped acrylic part for y-direction alignment
- oil coated (small friction) → avoid glue adhering




- an acrylic part for x & z direction alignment





◇ Assembling


- Optical cement
 - used for coupling between PCBs and scintillation bars
 - aligned with the jig



the smallest counter




- Aluminized mylar
 - used as a reflector



Counter Test

- Time resolution of assembled counters




the smallest counter was used as a reference counter

$\rightarrow \sigma \sim 65 \text{ psec}$

- $[T(\text{ch1}) + T(\text{ch2})] / 2 =: T_{\text{mean}}$
- $[T(\text{ch3}) + T(\text{ch4})] / 2 =: T_{\text{ref}}$



$T_{\text{mean}} - T_{\text{ref}}$ distribution



	time resolution
15 cm	71 psec
19 cm	104 psec

\rightarrow good enough for RDC

a frame for plastic scintillators



◇ Summary

- 16 ~ 28% sensitivity improvement will be achieved by introducing RDC.
- Downstream RDC has already been approved.
- Plastic scintillator bars have arrived and showed good time resolution.
- Counter production has started.
- Assembled counters showed good performance.

◇ Prospect

- Finish counter production and their test
- Combine plastic scintillator part with LYSO part, and test them as overall downstream RDC.
 - Downstream RDC will be prepared by the end of this year.
- Study on Upstream RDC with beam will be started this autumn.

END

BACKUP

◇ Plastic scintillator

Physical Constants of SGC Plastic Scintillators





Scintillator	Light Output % Anthracene ¹	Wavelength of Maximum Emission, nm	Decay Constant, Main Component, ns	Bulk Light Attenuation Length, cm	Refractive Index	H:C Ratio	Loading Element % by weight	Density	Softening Point °C
BC-400	65	423	2.4	250	1.58	1.103		1.032	70
BC-404	68	408	1.8	160	1.58	1.107		1.032	70
BC-408	64	425	2.1	380	1.58	1.104		1.032	70
BC-412	60	434	3.3	400	1.58	1.104		1.032	70
BC-416	38	434	4.0	400	1.58	1.110		1.032	70
BC-418	67	391	1.4	100	1.58	1.100		1.032	70
BC-420	64	391	1.5	110	1.58	1.102		1.032	70
BC-422	55	370	1.6	8	1.58	1.102		1.032	70
BC-422Q	11	370	0.7	<8	1.58	1.102	Benzophenone,0.5%*	1.032	70
BC-428	36	480	12.5	150	1.58	1.103		1.032	70
BC-430	45	580	16.8	NA	1.58	1.108		1.032	70
BC-440	60	434	3.3	400	1.58	1.104		1.032	99
BC-440M	60	434	3.3	380	1.58	1.104		1.039	100
BC-444	41	428	285	180	1.58	1.109		1.032	70
BC-452	32	424	2.1	150	1.58	1.134	Lead,5%	1.080	60
BC-454	48	425	2.2	120	1.58	1.169	Boron,5%	1.026	60
BC-480	**	425	—	400	1.58	1.100		1.032	70
BC-482A	QE=.86	494	12.0	300	1.58	1.110		1.032	70
BC-490	55	425	2.3	NA	1.58	1.107		1.032	70
BC-498	65	423	2.4	NA	1.58	1.103		1.032	70

¹ Anthracene light output = 40-50% of NaI(Tl)

* 0.1 to 5 weight % also available

** Ratio of Cerenkov light to scintillator light = 10:1

COUNTER PRODUCTION & TEST



◇ MPPCs for plastic scintillator

- MPPC (Multi-Pixel Photon Counter)
 - S13360-3050PE
 - pixel size: $50 \times 50 \mu\text{m}^2$
 - crosstalk suppression
 - 2 or 3 series connection
 - 60 MPPCs are used



trenches between pixels suppress crosstalk